

July 2024

Dear Customer,

IMPORTANT NOTICE: Glove Plus Ultra Purple – Ammonium Hydroxide 25%

As part of ongoing compliance with the requirements of EU Regulation 2016/425 on Personal Protective Equipment (PPE), product on the market is tested by our notified Body SATRA Limited.

During this routine testing, SATRA have advised us of an error on our packaging regarding the Permeation Resistance of Glove Plus Ultra Purple gloves (Product Code UNP3) in contact with:

- **O - Ammonium Hydroxide 25%.**

Our packaging states the Permeation Resistance (EN ISO 374-1:2016+A1:2018/Type B) is Level 5. This is incorrect.

The correct level of Permeation Resistance (EN ISO 374-1:2016+A1:2018/Type B) for Glove Plus Ultra Purple in contact with Ammonium Hydroxide 25% is **Level 2**.

Our packaging has been updated and will display the correct value from December 2024 onwards. The product data sheet has been revised and is available for download now in the 'Document Downloads' section of this web page. The latest version of the product data sheet is: PDS13-Eng-Jul 24.


The reduction in chemical permeation performance of Glove Plus Ultra when in contact with Ammonium Hydroxide 25% does not affect the 'type' categorisation of the product – this remains unchanged at **Type B**.

The *Determination of resistance to degradation by chemical perforation test (% degradation EN ISO 374-4:2019)* **16.1%** stated on our packaging and supporting documentation is correct.

Further information regarding the performance characteristics of Ultra Purple is available in the 'Document Downloads' section of the product page.

Should you require any further information, please do not hesitate to email: hello@barberhealthcare.com or call our team on +44 (0) 1969 624 846 (lines open Monday – Friday 8:00am to 5:00pm GMT).

Yours faithfully,



James Barber
Managing Director
Barber Healthcare Limited